

Neutron Cross Section Measurements at ORELA

K. H. Guber¹, L. C. Leal¹, J.A. Harvey²,
R. O. Sayer¹, H. Derrien¹, D. A. Wiarda¹
and P. E. Koehler²

¹Nuclear Science and Technology Division

²Physics Division

Oak Ridge Electron Linear Accelerator



- **High flux** (10^{14} n/sec) => gram-sized, affordable samples
- **Excellent resolution** ($\Delta t=2-30$ ns) => good S/N facilitates better evaluations
- **"White" neutron spectrum** from $E_n \sim 0.01$ eV - 80 MeV => reduces systematic uncertainties
- **Measurement systems and backgrounds well understood** => very accurate data
- **Simultaneous measurements** => (n,γ) , (n,α) , (n,n') , (n,f) , and σ_{total} experiments at the same time on different beam lines
- **Measurements on over 180 Isotopes:** ORELA measurements have contributed to ~80% of U.S. Evaluated Nuclear Data File (ENDF/B) evaluations




Neutron Cross-Section Measurements in Resolved Resonance Region

- Flight path length, **the longer the better**
- Pulse width of the neutron burst, **the shorter the better**:
 - Typically fixed with spallation sources (tens of ns to hundreds of ns)
 - Linac sources can vary pulse width (1ns up to tens of ns)
- Source moderation distance
 - Uncertainty of the creation location of the neutron inside the target/moderator configuration impacts resolution
 - ORELA has relatively small neutron production target
 - Spallation sources are usually optimized for thermal neutron flux, this requires large moderators
 - The moderation effect places tails on resonances due to delayed neutron emission
 - Hinders resolution of closely spaced resonances
 - Produces background in unresolved energy range which cannot be corrected. This effect is on the order of 16% at 20 keV (Coceva et al. 2002) for n_TOF and can not be estimated quantitatively

The Influence of Resolution



Doppler - Resolution Broadening @20 m for ^{241}Pu




Modified capture data measurement system has significantly less structure

Sample changer replaced


Aluminum guide replaced by C-fiber tube

Detector housing removed





Much of the Old Neutron Data (on Which Current Evaluations Are Based) Are Seriously Incorrect

- Some problems with the old data:
 - Underestimated neutron sensitivity correction
 - Low-energy cut off of 3 keV
 - No high energy (>100 keV) data
 - Incorrect weighting function
 - Poor resolution
 - Poorly characterized samples, i.e. water in the sample



Ex: Large neutron sensitivity of older measurements led to many erroneously-large resonance areas in current evaluations.




2D: Counts as a function of pulse-height and TOF

New Set-up for Transmission Measurements

- Time and pulse-height signals are digitized in two channels of an Acqiris digitizer
- Allows to take 2D spectra: Counts as a function of TOF and pulse-height
- Fewer electronic components, thus simpler and more reliable
- Unlimited stops per start
- Easy upgrade path to allow the use of ^6Li -glass and NE110 simultaneously in the same transmission measurement. This would allow to cover the whole energy range with the most efficient detector.

The Importance of Total Cross Section Data

- More complete resonance parameter data will help improve nuclear statistical model.
- Is indispensable for obtaining the most accurate (n,γ) reaction rates. **See resonances not visible in (n,γ) data. Improved self-shielding and multiple scattering corrections.**
- Lack of good total cross section data can lead to serious errors in these corrections and hence in the cross sections.
- **Ex: ^{116}Sn Use of incorrect neutron widths led to incorrect low-energy cross sections (Wisshak et al.).**




Elastic Scattering Measurement Capabilities

- Scattering Chamber at 200m with six neutron detectors
 - 4.32 cm diameter, 7.62 cm length NE110 cylinder
 - 8850 photomultiplier tubes connected at each end of NE110 cylinder
- Position detectors at various lab angles from target sample located in center of scattering chamber





^{58}Ni Measurement and SAMMY Analysis





Recent ORELA Measurements

- Use of a 99% enriched ^{41}KCl sample in the neutron beam of ORELA at flight path #7 in the 40 meter flight station
- Beam specifications: 525 Hz, 8 nsec pulse width and 4kW
- New electron gun design performing well
- Able to start/stop accelerator as desired—vacuum holding
- Further improvements in beam power anticipated as the system is fine-tuned



- Restarted ORELA April 17, 2006
- ~100 hours beam time in last couple months





Differential distribution of resonance spacing (Wigner distribution) for ^{238}U



Distributions Of The Reduced Neutron Widths (Porter-Thomas distribution) for ^{238}U




Neutron Capture and Total Cross Section Experiments at a White Neutron Source: Typical Results



Only small part (0.06%) of energy range of data shown.

Separate sample-out and background measurements needed.

Existing Experiments at ORELA



- 11 Flight paths
- Flight Stations:
 - 8-18, 20, 35, 40, 85, 150, and 200 m